Compuex®

Fit 5.0

SP 6.0

SP 8.0

Quick Start Guide

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Beknopte handleiding

Guia rápido

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
Error

Erreur • Fehler
Error • Errore
Fout • Erro
Error • Ошибка

**Contact the
Complex-approved customer
service center**

Contacter le service clientèle agréé par Complex
Wenden Sie sich an ein von Complex zugelassenes Service-Zentrum.
Contattare il centro di assistenza autorizzato Complex
Pongase en contacto con el centro de servicio al cliente aprobado por Complex
Neem contact op met het door Complex erkende klantenservicecentrum
Contacte o centro de assistência ao cliente aprovado pela Complex
обратитесь в авторизованный ремонтный центр Complex.

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


Well done! You may now also try other programmes

- Bravo !** Vous pouvez maintenant essayer d'autres programmes
- Gratuliere!** Sie können jetzt andere Programme ausprobieren
- Bien hecho!** Ahora debería probar otros programas
- Ben fatto!** Ora puoi provare altri programmi
- Goed gedaan!** U kunt nu ook andere programma 's proberen
- Correcto!** Agora pode experimentar outro programa
- Bra gjort!** Nu kan du även prova andra program
- Отлично!** Теперь вы можете попробовать другую программу

1. Position of electrodes

For optimum results, use the recommended electrode positions.

A different electrode must be connected to each extremity of the cable.

N.B.: It is possible and normal that large dual-connection electrodes are used in such a way that a connection is not connected.

It is important to select appropriately sized electrodes (small or large) and to position them correctly on the muscle group to be stimulated in order to ensure effective treatment. It is therefore appropriate to always use the electrode size shown on the illustration. If necessary, try to find the best position possible by moving the electrode on the muscle until the point that produces the best contraction or greatest comfort is found.

For TENS type programmes, the general rule is to cover the painful area with electrodes.

N.B.: Approximate positioning of electrodes makes the session less effective, but is not dangerous.

2. Body position during stimulation

This position varies according to the muscle group to be stimulated and the programme that you are using.

For programmes inducing powerful muscular contractions (tetanic contractions), it is recommended stimulating the muscle isometrically. Therefore, it is necessary to hold the extremities of limbs stimulated. This position opposes maximum resistance to movement and avoids any shrinking of muscle during contraction, which would risk causing cramps. For example, when stimulating quadriceps, you should be in a sitting position with ankles held to avoid extending the knees.

It is also possible to perform dynamic effort associated with your Compex i.e. perform a voluntary dynamic movement at the same time as the stimulation; for example, perform squats during stimulation of the quadriceps. But this way of working is reserved for expert users who are accustomed to the technique of electrostimulation who have already performed complete stimulation cycles (minimum 4 to 6 weeks at a rate of 3 sessions /week). In addition, the introduction of dynamic effort must be very progressive; for example 10 contractions the first week (the remainder of the session in isometric), then 15 the following week, etc.

In this case, it is also important to perform safe and simple dynamic movements because the addition of electrostimulation to dynamic effort significantly increases the workload imposed on the muscle. For example, avoid performing dynamic movements with impact (e.g. jumping from a bench) or adding heavy loads such as weights during exercise.

For the other types of programmes which do not induce powerful muscular contractions but only muscular twitches or tingling, position yourself as shown on the illustration but always in a comfortable position.

3. Adjusting stimulation intensities

On a stimulated muscle, the number of fibres used depends on the stimulation intensity.

For programmes inducing powerful muscular contractions (tetanic contractions), the maximum stimulation intensities must be used (up to 999) always at the limit of what you can tolerate, in order to use a maximum number of fibres. The progress of a muscle will be greater the higher number of fibres that are worked. During a session it is therefore important to try to increase the intensity every 3 to 5 contractions.

For other types of programmes like recovery, massage, capillarisation or even muscle pain which only induce muscle twitches, gradually increase the stimulation intensities until visible muscle twitches are obtained.

For TENS, Epicondylitis and Tendinitis programmes, gradually increase the intensities until tingling sensation are felt under the electrodes.

4. Progression of stimulation cycles

The stimulation cycles are for people who are already accustomed to electrostimulation and want to perform several training cycles.

The cycle logic refers to the workload performed by electrostimulation. And just like a normal workout, one has to start with an amount of effort then increase it over the course of the cycles. Thus, it is recommended starting with the 1st cycle and going on to the next level when the cycle is finished, normally after 4 to 6 weeks of stimulation based on 3 sessions per week. It is also important to have reached significant stimulation intensities in sessions before going on to another cycle.

At the end of a cycle you can either start a new cycle or do maintenance training based on one session per week.

N.B.: To be accessible, the Cycles function must be enabled in the device's settings.

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